

CENTOGENE
THE RARE DISEASE COMPANY



CentoCancer[®]

Strive for the
complete information

Centocancer® – Panel Composition and Methodology

Centocancer® includes the following 70 most relevant cancer associated genes:

<i>ABRAXAS1</i>	<i>BRIP1</i>	<i>FLCN</i>	<i>MLH3</i>	<i>PMS1</i>	<i>RAD51D</i>	<i>SMARCA4</i>
<i>APC</i>	<i>CDH1</i>	<i>GALNT12</i>	<i>MRE11</i>	<i>PMS2</i>	<i>RECQL</i>	<i>STK11</i>
<i>ATM</i>	<i>CDK4</i>	<i>HNF1B</i>	<i>MSH2</i>	<i>POLD1</i>	<i>RET</i>	<i>TGFBR2</i>
<i>AXIN2</i>	<i>CDKN2A</i>	<i>HOXB13</i>	<i>MSH3</i>	<i>POLE</i>	<i>RNF43</i>	<i>TP53</i>
<i>BAP1</i>	<i>CHEK2</i>	<i>KIT</i>	<i>MSH6</i>	<i>POT1</i>	<i>SDHA</i>	<i>TSC1</i>
<i>BARD1</i>	<i>DICER1</i>	<i>MC1R</i>	<i>MUTYH</i>	<i>PRSS1</i>	<i>SDHAF2</i>	<i>TSC2</i>
<i>BLM</i>	<i>DIS3L2</i>	<i>MEN1</i>	<i>NBN</i>	<i>PTCH1</i>	<i>SDHB</i>	<i>VHL</i>
<i>BMPR1A</i>	<i>EPCAM</i>	<i>MET</i>	<i>NF1</i>	<i>PTEN</i>	<i>SDHC</i>	<i>WT1</i>
<i>BRCA1</i>	<i>FANCC</i>	<i>MITF</i>	<i>NTHL1</i>	<i>RAD50</i>	<i>SDHD</i>	<i>XRCC2</i>
<i>BRCA2</i>	<i>FH</i>	<i>MLH1</i>	<i>PALB2</i>	<i>RAD51C</i>	<i>SMAD4</i>	<i>XRCC3</i>

Key panel facts

- Next-generation sequencing (NGS) of all 70 genes in the panel, including all coding regions and +/-10bp exon/intron boundaries
- Coverage: ≥ 99.5% of target bases covered at > 20x
- NGS-based CNV (copy number variant) analysis for all genes
- Low quality single nucleotide variants (SNVs) and all relevant deletion/insertion variants are confirmed by Sanger sequencing or MLPA/qPCR prior to reporting
- All relevant deep intronic variants described in the current version of HGMD® and our rare disease-centric Bio/Databank are included
- Turnaround Time: 15 business days
- Required Material: ≥ 1 µg DNA or ≥ 1 ml EDTA blood or ≥ 1 CentoCard®

Centocancer – Our Comprehensive Oncogenetics Panel for Hereditary Mutations

Hereditary pathogenic variants confer an increased risk of developing cancers during an individual's lifetime. Early identification of pathogenic variants in genes which have a predisposition to cancer is a fundamental first step in the diagnosis, management and treatment of individuals and families with hereditary cancer syndromes.

Panel Composition

Centocancer offers complete answers to help you choose the best possible therapeutic approach for your patients. Each gene in Centocancer has been carefully selected based on its risk potential in the development of one or more of the following cancers:



Breast



Colorectal



Thyroid



Pancreatic



Renal



Ovarian



Gastric



Endometrial



Melanoma




Prostate

Who should consider Centocancer for genetic testing?

Centocancer is appropriate for:

- 1 Individuals with a positive personal history of early-onset cancer, rare cancer, bilateral cancer, or multiple primary cancers
- 2 Unaffected individuals with a positive family history of multiple generations of cancers, rare cancers, or early-onset cancers
- 3 Individuals in whom the suspected genetic diagnoses for a suspected familial cancer risk are not covered by a single targeted panel, or if a targeted panel testing was previously negative



Hereditary Cancer and/or Susceptibility

SELECTION OF GENETIC TEST/PANEL ACCORDING TO FAMILY HISTORY AND CLINICAL DATA

COMPLEX FAMILY HISTORY, VARIABILITY OF CANCERS AND ABSENCE OF KNOWN GENETIC CAUSE IN THE FAMILY

BRCA1, BRCA2 panel	<i>BRCA1, BRCA2</i>
CentoBreast®	<i>ABRAXAS1, ATM, BARD1, BRCA1, BRCA2, BRIP1, CDH1, CHEK2, DICER1, EPCAM, FANCC, MEN1, MLH1, MRE11, MSH2, MSH6, MUTYH, NBN, PALB2, PMS1, PMS2, PTEN, RAD50, RAD51C, RAD51D, RECQL, SMARCA4, STK11, TP53, XRCC2</i>
CentoColon	<i>APC, ATM, AXIN2, BLM, BMPR1A, BRCA1, BRCA2, CDH1, CDKN2A, CHEK2, EPCAM, FLCN, GALNT12, MLH1, MLH3, MSH2, MSH3, MSH6, MUTYH, NBN, NTHL1, PALB2, PMS2, POLD1, POLE, PTEN, PRSS1, RNF43, SMAD4, STK11, TGFBR2, TP53, VHL</i>

CentoCancer panel	<i>ABRAXAS1, APC, ATM, AXIN2, BAP1, BARD1, BLM, BMPR1A, BRCA1, BRCA2, BRIP1, CDH1, CDK4, CDKN2A, CHEK2, DICER1, DIS3L2, EPCAM, FANCC, FH, FLCN, GALNT12, HNF1B, HOXB13, KIT, MC1R, MEN1, MET, MTF, MLH1, MLH3, MRE11, MSH2, MSH3, MSH6, MUTYH, NBN, NFI, NTHL1, PALB2, PMS1, PMS2, POLD1, POLE, POT1, PRSS1, PTCH1, PTEN, RAD50, RAD51C, RAD51D, RECQL, RET, RNF43, SDHA, SDHAF2, SDHB, SDHC, SDHD, SMAD4, SMARCA4, STK11, TGFBR2, TP53, TSC1, TSC2, VHL, WT1, XRCC2, XRCC3</i>
--------------------------	---

Identification of specific cancer-causing pathogenic variant

No pathogenic variants identified

Identification of specific cancer-causing pathogenic variant

No pathogenic variant identified

WES analysis on a research basis

Research reporting

Genetic counseling, genetic testing of all family members with consent

Genetic counseling, genetic testing of all family members with consent

Some Common Cancer Predisposition Syndromes Covered by CentoCancer

Syndromes

HEREDITARY BREAST/OVARIAN CANCER

BRCA1, BRCA2

LI-FRAUMENI SYNDROME

TP53

COWDEN SYNDROME

PTEN

HNPCC (LYNCH SYNDROME)

MLH1, MSH2, MSH6, PMS1, PMS2

FAMILIAL ADENOMATOUS POLYPOSIS

APC

VON HIPPEL-LINDAU

VHL

MULTIPLE ENDOCRINE NEOPLASIA

MEN1, RET

Associated Cancers

- > Breast, ovarian, prostate, pancreatic, melanoma
- > Breast, sarcomas, adrenocortical carcinoma, leukemia, brain tumors
- > Breast, thyroid, benign lesions of skin, hamartoma, renal cell carcinoma, uterine
- > Colorectal endometrial, ovarian, small bowel, stomach, pancreas, ureter, renal pelvis
- > Polyposis, colorectal, thyroid, gastric, periampullary carcinoma, hepatoblastoma
- > Renal cell carcinoma, retinal angioma, cerebellar hemangioblastoma, pheochromocytoma, pancreatic cysts, islet cell tumor
- > Parathyroid tumors, pancreatic tumors, pituitary tumors, medullary thyroid cancer, pheochromocytoma, neuromas

Please visit our website for more information:
www.centogene.com

CONTACT DETAILS

CENTOGENE GmbH

Am Strande 7
18055 Rostock
Germany

CENTOGENE GmbH is a subsidiary of CENTOGENE N.V.

✉ customer.support@centogene.com

☎ +49 (0) 381 80 113 - 416

📄 +49 (0) 381 80 113 - 401

Rostock - CLIA #99D2049715



Cambridge - CLIA #22D2154474

