

# Solid Tumor Panel

Providing Knowledge to Battle Cancer

PRODUCT SHEET

## Solid Tumor Panel

Somatic mutations are acquired genetic modifications unique to an individual that can affect any type of cell in the body, except for the germ cells. Somatic mutations often also consist of known 'driver mutations,' which sometimes can be targeted by highly specific and efficient drugs. The knowledge of those targetable mutations in tumor cells can therefore identify additional treatment options.

CENTOGENE's Solid Tumor Panel with a total of 149 genes provides full sequencing of 106 selected cancer-associated genes as well as the hotspot analysis of relevant cancer regions in 43 genes. It detects over 5,000 validated oncogenic variants and includes the latest evidence-based variants associated with treatment decisions in solid tumors. The panel covers as more than 25 genetic variants with directed therapies which follow the FDA and NCCN guidelines as well as many that are currently being tested in clinical trials. Furthermore, somatic variants with an impact on prognosis of the individual tumor or on the efficacy of standard anti-tumor therapy are captured and reported.

### The CENTOGENE Advantage



Curated to **provide the most valuable information** for diagnostic decisions, prognosis, and therapeutic approaches



Powered by CENTOGENE's Biodatabank, **the world's largest real-world data repository** for rare and neurodegenerative diseases



The most **up-to-date panel gene content** with the latest medical and in-house findings



**High-quality analysis for precise clinical interpretation** using advanced bioinformatics and artificial intelligence-powered tools

## Genes Included (149)

Coverage of full coding regions for:

*ABL1, AKT1, AKT2, AKT3, APC, AR, ARID1A, ASXL1, ATM, ATR, ATRX, BAP1, BRAF, BRCA1, BRCA2, CDH1, CDK12, CDK4, CDKN1B, CDKN2A, CDKN2B, CHEK1, CHEK2, CREBBP, CSF1R, CTNNA1, DDR2, EGFR, ERBB2, ERBB3, ERBB4, EZH2, FANCA, FANCD2, FANCD2, FANCI, FBXW7, FGFR1, FGFR2, FGFR3, FGFR4, GNA11, GNAQ, GNAS, HNF1A, HRAS, IDH1, IDH2, KDR, KEAP1, KIT, KMT2A, KMT2C, KMT2D, KRAS, MAO2K1, MAP2K2, MEN1, MET, MLH1, MPL, MRE11, MSH2, MSH6, MTOR, NBN, NF1, NF2, NFE2L2, NOTCH1, NOTCH2, NOTCH3, NRAS, NTRK3, PALB2, PDGFRA, PIK3CA, PIK3R1, PMS2, POLE, PTCH1, PTEN, PTPN11, RAD50, RAD51, RAD51B, RAD51C, RAD51D, RB1, RBM10, RET, RIT1, RNF43, SETD2, SLX4, SMAD4, SMARCA4, SMARCB1, SMO, SPOP, SRC, STK11, TP53, TSC1, TSC2, TSHR, VHL* with +/- 2bp flanking intronic regions

Hotspot analysis targeting relevant cancer-associated regions for the following genes:

*ALK, ARAF, AXL, BTK, CBL, CCND1, CDK6, ERCC2, ESR1, FLT3, FOXL2, GATA2, H3-3A, H3C2, JAK1, JAK2, JAK3, KNSTRN, MAGOH, MAP2K4, MAPK1, MAX, MDM4, MED12, MYC, MYCN, MYD88, NTRK1, NTRK2, PDGFRB, PIK3CB, PPP2R1A, RAC1, RAF1, RHEB, RHOA, ROS1, SF3B1, STAT3, TERT, TOP1, U2AF1, XPO1*

## Key Features and Performance

<b>COVERAGE</b>	≥ 97.0% targeted regions covered at ≥ 200x
<b>VARIANT TYPES</b>	<ul style="list-style-type: none"> <li>• Sensitivity SNVs and InDels (≤ 50bp)* &gt; 97.2%</li> <li>• Accuracy of &gt; 96.2%</li> <li>• Specificity of ≥ 99.9% guaranteed for all reported variants. Variants with low quality and/or unclear zygosity are confirmed by orthogonal methods**</li> </ul>
<b>REPORTING</b>	Pathogenic and likely pathogenic variants are reported following ACMG classification guidelines recommendations. Additionally, these variants are reported according to their actionability into Tier 1 (strong clinical significance) or Tier 2 (potential clinical significance), following the standards and guidelines for the Interpretation and reporting of sequence variants in cancer.***
<b>REQUESTED MATERIAL</b>	FFPE tissue (block or sections)**** or fresh tumor tissue
<b>TAT</b>	10 business days

SNVs: single nucleotide variants; InDels: small insertions/deletions

\* SNV detection down to 5% allele frequency

\*\* Variants with low quality and / or unclear zygosity are confirmed by orthogonal methods, i.e.: SNVs and InDels by Sanger sequencing

\*\*\* Li et al. 2017, PMID: 27993330

\*\*\*\* Formalin-Fixed Paraffin Embedded (FFPE) Block or ≥10 sections of 10µm thickness with control slide

## Common Directed Therapies Covered by Solid Tumor Panel

BIOMARKER*	MECHANISM	EVIDENCE LEVEL/ GUIDELINES
ABL1 (T315I)	BCR-ABL inhibitor 3rd gen&Pan-TK inhibitor	FDA
BRAF (V600E)	BRAF inhibitor	FDA
BRAF (V600E,V600K)	BRAF inhibitor;MEK inhibitor	FDA
BRAF (V600E,V600K)	MEK inhibitor	FDA
BRAF (V600E)	BRAF inhibitor	FDA
BRAF (V600E,V600K)	BRAF inhibitor;MEK inhibitor	FDA
BRCA1 oncogenic mutations (complete gene covered)	PARP inhibitor	FDA
BRCA2 oncogenic mutations (complete gene covered)	PARP inhibitor	FDA
EGFR exon 19 deletions	ERBB2 inhibitor&EGFR inhibitor 2nd gen	FDA
EGFR (L858R,L861Q,G719A,G719S,G719C,G719D,L747S,S768I,L861P,L861Q,L861R)	ERBB2 inhibitor&EGFR inhibitor 2nd gen	FDA
EGFR exon 19 deletions	EGFR inhibitor 1st gen	FDA
EGFR (L858R)	EGFR inhibitor 1st gen	FDA
EGFR exon 19 deletions	EGFR inhibitor 1st gen	FDA
EGFR (L858R,L861Q,G719A,G719S,G719C,G719D,L747S,S768I,L861P,L861Q,L861R)	EGFR inhibitor 1st gen	FDA
EGFR (T790M)	EGFR inhibitor 3rd gen	FDA
JAK2 (V617F)	JAK inhibitor	FDA
KIT mutation in exon 9,11,13,14 or 17	BCR-ABL inhibitor 1st gen&KIT inhibitor	FDA
KIT mutation in exon 9,11,13,14 or 17	Pan-kinase inhibitor	FDA
KIT mutation in exon 9,11,13,14 or 17	Pan-TK inhibitor	FDA
PTCH1 oncogenic mutations (complete gene covered)	SHH inhibitor	FDA
RET (618,620,634,768,791,891,918,C634W,M918T)	Pan-TK inhibitor	FDA
TSC1 oncogenic mutations (complete gene covered)	MTOR inhibitor	FDA
TSC1 oncogenic mutations (complete gene covered)	MTOR inhibitor	FDA
TSC2 oncogenic mutations (complete gene covered)	MTOR inhibitor	FDA
ABL1 (T315A,F317L,F317V,F317I,F317C,F317I,Y253H,E255K,E255V,F359V,F359C,F359I)	BCR-ABL inhibitor 3rd gen	NCCN
ABL1 (F359V,F359C,F359I,Y253H,E255K,E255V)	BCR-ABL inhibitor 2nd gen	NCCN
ABL1 (T315A,F317L,F317V,F317I,F317C,V299L)	BCR-ABL inhibitor 2nd gen	NCCN
ABL1 (T315I)	BCR-ABL inhibitor 3rd gen&Pan-TK inhibitor	NCCN
BRAF (V600D,V600K,V600M,V600G,V600R)	BRAF inhibitor	NCCN
BRAF (V600E,V600D,V600K,V600M,V600G,V600R)	BRAF inhibitor	NCCN
EGFR (L861Q,G719A,G719S,G719C,G719D,L747S,S768I,L861P,L861Q,L861R)	ERBB2 inhibitor&EGFR inhibitor 2nd gen	NCCN
EGFR (L858R,L861,G719,S768I)	EGFR inhibitor 1st gen	NCCN
EGFR (L861Q,G719A,G719S,G719C,G719D,L747S,S768I,L861P,L861Q,L861R)	EGFR inhibitor 1st gen	NCCN
KIT mutation in exon 9,11,13,14 or 17	BCR-ABL inhibitor 1st gen&KIT inhibitor	NCCN
PDGFRA inframe deletion (I843)	BCR-ABL inhibitor 1st gen&KIT inhibitor	NCCN
PDGFRA (552-596,631-668,814-854)	BCR-ABL inhibitor 1st gen&KIT inhibitor	NCCN
PDGFRA (552-596,631-668,814-854)	Pan-kinase inhibitor	NCCN
PDGFRA (552-596,631-668,814-854)	Pan-TK inhibitor	NCCN

\*based on the CGI Cancer Biomarkers database (2018/01/17 data)