


**CENTOGENE**  
THE RARE DISEASE COMPANY



**Centoswab™**

**YOUR SAMPLE COLLECTION SOLUTION  
TO DETECT SARS-COV-2 RNA**

Testing of SARS-CoV-2 RNA in the midst of a  
global shortage

# Centoswab™

FACILITATING FAST AND RELIABLE SARS-COV-2 RNA TESTING

## Key Benefits




- Extensively validated with patient and control samples for the detection of SARS-CoV-2 RNA
- INSTAND e.V. quality assurance for SARS-CoV-2 testing with inactivated positive and negative proband material
- Ethylene-oxide treated, ready to use swabs and collection tubes
- Medical device class I – CE labeled
- Easy logistics – can be sent without stabilization buffer
- Made in Germany

## Sample Stability - E gene SARS-CoV-2 RNA



- Collected sample RNA on Centoswab™ is stable for 72 h after collection<sup>1</sup>
- Centoswab™ has an excellent sensitivity of < 2 virus particles as demonstrated with validated SARS-CoV-2 RT-PCR assay



## Centoswab™ Kit Includes



- Oropharyngeal Swab
- Collection Tube

## ➤ Contact Details

Phone: +49 (0)381 80 113-700

Email: [centoswab@centogene.com](mailto:centoswab@centogene.com)

[www.centogene.com](http://www.centogene.com)

CLIA #99D2049715



<sup>1</sup> Sample must be stored and shipped at 2 – 8°C after collection if not processed immediately

<sup>2</sup> For further test details please visit our homepage [www.centogene.com/corona](http://www.centogene.com/corona)